

NOW AVAILABLE!

Dyanavel[®]XR[®]
 (amphetamine) extended-release
 oral suspension 2.5 mg/mL

**New
60 mL
Bottle**



SMALLER SIZE for convenient pharmacy stocking

- ✓ For the treatment of ADHD in patients 6 years and older
- ✓ **Bubble-gum flavored** oral suspension

How Supplied

DYANAVEL[®] XR (amphetamine) extended-release oral suspension: The concentration is 2.5 mg/mL amphetamine base equivalents and is supplied as light beige to tan viscous suspension with bubblegum flavor in bottles of 60 mL (**NDC 24478-102-02**) and 464 mL (**NDC 24478-102-01**).

The product is provided in a carton.


- 60 mL carton contains two oral dispensers and two bottle adapters.
- 464 mL carton also contains four oral dispensers and four bottle adapters.

Wholesaler		Product Description
Cardinal	CIN# 6108096	DYANAVEL XR (amphetamine) Extended-Release Oral Suspension 2.5 mg/mL
Cencora	1037665	DYANAVEL XR (amphetamine) Extended-Release Oral Suspension 2.5 mg/mL
McKesson	307665	DYANAVEL XR (amphetamine) Extended-Release Oral Suspension 2.5 mg/mL

Dosage Dispensing Instructions

Daily dosage	2.5 mg	5 mg	7.5 mg	10 mg	12.5 mg	15 mg	17.5 mg	20 mg
DYANAVEL XR volume	1 mL	2 mL	3 mL	4 mL	5 mL	6 mL	7 mL	8 mL
Dispense (30 day supply)	30 mL	60 mL bottle	60 mL bottle + 30 mL	2x 60 mL bottle	2x 60 mL bottle + 30 mL	3x 60 mL bottle	3x 60 mL bottle + 30 mL	4x 60 mL bottle

Please see accompanying Full Prescribing Information, including **BOXED WARNING** regarding **ABUSE, MISUSE, and ADDICTION**, or scan the QR code on the back.



Patient instructions

- Use with the oral dosing dispenser provided by the pharmacist
- Ensure that the bottle adapter has been firmly inserted into the bottle by the pharmacist. Once the bottle adapter has been inserted into the bottle, it should not be removed
- Shake the bottle of DYANAVEL XR before each dose
- Measure the appropriate dose as prescribed
- Using the filled oral dosing dispenser, dispense DYANAVEL XR oral suspension directly into mouth
- Replace bottle cap and store bottle as directed
- Wash oral dosing dispenser after each use



Scan for dosing calculator or visit
dyanavelxrhcp.com/dyanavel-xr-liquid/dosing

Storage and Handling

Dispense in a tight, light-resistant container with child-resistant closure.

Store at 20° to 25°C (68° to 77°F); excursions permitted from 15° to 30°C (59° to 86°F) [see USP Controlled Room Temperature].

DYANAVEL XR extended-release oral suspension:

The pharmacist should insert the bottle adapter firmly into the neck of the bottle and provide the oral dosing dispenser to the patient when dispensing this product.



Please see accompanying Full Prescribing Information, including BOXED WARNING regarding ABUSE, MISUSE, and ADDICTION, or scan the QR code.